

MVAK

DE

GB

IT

CE

Inhalt / Content / Contenuto

Seite / Page / Pagina

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Zubehör

Differenzdruck-Manometer
"Econometer"

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

Functional elements

- A Econometer ¹⁾
Differential pressure manometer
- B Sealing plate
- C Flat seal
- D Upper housing bowl
- E Housing O-ring
- F Filter element(s)
- G Lower housing bowl
- H Coupling nut
- I Ball valve and drain flask ²⁾
- J Tie rod and knurled nut
- K Sight glass and seal
- ^{1), 2)} according to configuration, see separate documentation

Attention!
Biological risk

For your safety

The relevant safety at work and accident prevention regulations, plus the operating instructions, shall apply for operation of the pressure vessel. The pressure vessel has been constructed in accordance with the generally recognized rules of engineering. It complies with the requirements of directive 97/23/EC concerning pressure vessels.

The relevant applicable national regulations in force at the place of installation concerning the operation and routine testing of pressure vessels must be complied with.

You as operator / user of the unit should make yourself familiar with the function, installation and start-up of the unit through these operating instructions.

It is essential that you follow these safety notes and this information in order to ensure trouble-free operation of the unit.

All the safety information is always intended to ensure your personal safety!

- The max. working pressure and the max. permissible working temperature of the pressure vessel are detailed on the type plate.
- The permissible working temperatures for add-on parts and filter elements are given under Technical data in these instructions.
- Ensure that the permitted operational temperatures are complied with, regardless of the ambient temperatures prevailing at the place of installation.
- It is necessary to ensure that the unit is equipped with the corresponding safety and test devices to prevent the permissible operating parameters from being exceeded.
- The pressure vessel must be at a safe distance of min. 5 m to prevent heating up in the event of a fire.

- The pressure vessel has been designed for a primarily static pressure loading with a maximum of 1000 cycles to and from the full load. Rapid changes of load with more than 10 % of the max. working pressure are not allowed.
- Ensure that the pressure vessel is not subjected to vibrations that could cause fatigue fractures.
- The pressure vessel is not to be subjected to stresses arising from traffic, wind and earthquakes.
- The medium used may not have any corrosive components that could attack the materials of the pressure vessel in a way that is not permitted.
- All Installation and maintenance work on the pressure vessel may only be carried out by trained and experienced specialists.
- It is forbidden to carry out any kind of work on the pressure vessel and piping, this covering welding and constructional changes, etc. Breaking this rule means extreme danger for you and your colleagues.
- Attention! If the pressure vessel is operated at temperatures over 60°C, suitable protection to prevent contact must be provided.
- A pressure gauge that shows the operational pressure must be installed in the unit.
- Depressurize the system before carrying out any work on the pressure vessel.
- Clean the piping before carrying out the installation work.
- The unit must be installed vertically in the piping.
- Ensure that the pressure vessel is installed without any stresses.
- Disconnect the power supply when carrying out electrical work.

Appropriate use

The equipment may only be used for its intended purpose. The equipment has been built exclusively:

- for operating media of group 2 as per Pressure Equipment Directive 97/23/EC.
- for the complete removal of bacteria, for the protection against contamination for vacuum-pumps and atmosphere.

Any other form of use or one going beyond this shall be considered as inappropriate. We shall have no liability whatsoever for any damage incurred as a result!

Notes on starting up

Before initial commissioning:

Take note of the flow direction!

The flow direction of the filter element is from outside to inside.

All the screwed connections of the pressure vessel must be done up to the required and max. permissible tightening torques for the screws and bolts.

Make a visual check! There must be no external damage visible.

Your unit can be equipped with various components. Note the permissible working pressures and permissible temperatures for the components. You can find information in the relevant operating instructions and also from the type plate of your unit or its components.

Make a check for leaks!

Initial commissioning:

Slowly apply pressure to the system by first opening the downstream valve (**Picture 11**).

Information concerning maintenance

Note!

In order to make the maintenance of the built-in filter easier and to avoid an interruption of the vacuum-supply, it is recommended to install two medical vacuum filters parallel.

In accordance with instructions wear protective clothes when maintenance work is being done, in order to avoid a biological contamination!

Before starting any maintenance work, ensure that the pressure vessel has been depressurized and has cooled down, and cannot be put back into operation during the maintenance work.

The filter elements must be changed at regular intervals. At the latest, once the permissible differential pressure has been reached!

A blocked filter element can lead to damage the vacuum-pump!

If the housing of your unit has been equipped with an Econometer, this will show optically (**Picture 12**) if the differential pressure has been exceeded. Otherwise, the following recommendation apply:

FFMV: Depending on the degree of dirtiness, change after one year at the latest.

The drain flask is to be controlled on a daily basis and emptied in regular intervals. Before renewed installation it must be cleaned and autoclaved.


The housing O-ring should also be changed at the same that the filter is changed (**Picture 18**).

Damaged components are to be replaced by new ones. If a marked degree of damage is found, the entire vessel is to be replaced.

The pressure vessel has been designed for a life of 10 years.

Carry out a check for leaks once the maintenance work has been finished!

Protection of the environment


Attention! In accordance with instructions a careful waste disposal of the used filter elements is imperative due to the basis of the biological contamination danger.

The packing material and the unit itself and its accessories are produced from recyclable materials.


Separating the remaining materials in an appropriate way helps in the recycling of materials.

Spare parts

Please always quote the type designation on your filter housing when making orders for spare parts.

Accessories

Please always quote the type designation on your filter housing when making orders for accessory parts.

Technical data

Max. working pressure PS: 2 bar a

Design temperature of the housing: - 25°C/+120°C

Permissible working temperature with filter type:

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


WORLDWIDE WARRANTY STATEMENT

Under proper use, installation, application, and maintenance in accordance with the manufacturer's written recommendations and specifications for a specific time period outlined below each ultrafilter product is warranted to be free from defects in materials and workmanship. Before installing the ultrafilter products, the customer has to ensure that his application/usage of the ultrafilter product is in accordance with the manufacturer's recommendations and specifications.

The manufacturer's warranty is a legal remedy and limited to the repair or replacement (at ultrafilter's option) of unaltered products returned to the manufacturer and proven to have such defect, provided the defects are reported within the time guideline listed below. Any replacement pads or dryers will extend to the remaining period of original warranty and not beyond.

Refrigerated Air Dryers and Chillers - 24 months from date of invoice.

Note: In addition to the above stated warranty, both air to air and air to refrigerant heat exchangers in ultrafilter dryers are warranted on a prorated basis from time of shipment for a period of five (5) years. This extended warranty is limited to these pads only and will be prorated in the following manner:

Year 1 - 100%, year 2 - 100 %, year 3 - 75% , year 4 - 50% , year 5 - 25% .

Regenerative Dryers - 24 months from date of invoice. Proof of installation must be received from claimant to verify installation date or commissioning by ultrafilter service.

All Other Commodity Products - filters, drains, economizers, ultraseps, and components, etc.: 12 months from date of final shipment to the end customer by factory or manufacturer's authorized distributor. 60 month material safety guarantee for housings from date of invoice.

ultrafilter original Spare parts- twelve (12) months ,pads only' warranty from date of purchase.

Products found to be damaged by exposure to acidic or otherwise corrosive environments will not be covered under normal or prorated warranty.

Maintenance and Adjustment - Adjustment to hot gas bypass valves, float drains, condenser cleaning, and other routine maintenance required to keep a dryer in good operating order are not considered a warranty item. This is the responsibility of the owner of the equipment, and labour or material for such will not be reimbursed under warranty. The usual maintenance and replacement type products are not covered by this warranty. Please consult the installation and operation manual for adjustment and maintenance requirements.

Before any warranty service work is started, it must first be authorized by ultrafilter. ultrafilter will not reimburse warranty work unless it has been preauthorized. Unauthorized service voids the warranty and any resulting charge or subsequent claim will not be paid.

What is And What Is Not Covered Under This Warranty?

NOT COVERED:

- All freight damage claims are not the responsibility of the manufacturer and are not covered under warranty, as all products are shipped F.O.B. shipper. Please direct freight claims to the shipper in question.
- More than one service technician on one warranty job, extra helpers deemed necessary (on-site personnel should assist as much as possible).
- Overtime hours, weekend labour rates.

- Airfreight of units or parts. If airfreight is required, the manufacturer will pay the cost for normal ground transportation and the customer will be liable for costs exceeding those charges.
- Removal or reinstallation of equipment or extra bypassing to accomplish repairs.
- Hot gas valve adjustments or other normal maintenance items.
- Damages due to misapplication, incorrect installation or use, poor maintenance, corrosive environments, or other items outside the scope of manufacturer defects. See operating and instruction books.
- Service work by unauthorized service.

COVERED:

- Repairs due to defects in materials.
- Repairs due to defects in workmanship.
- A maximum of 200 km round trip travel for one service technician at a maximum rate of 50 € / hour and 30 cent / per mile.
- Any accessories, pads and equipment supplied by ultrafilter, but not manufactured by ultrafilter shall carry whatever warranty the manufacturer has given to ultrafilter, provided it is possible for ultrafilter to pass on such warranty to the customer.
- Total value of claim not to exceed the value of the pad including labour and materials,

Specifications, limitations and recommended applications and uses for and of products may be established or amended by ultrafilter from time to time. Only those specifications, limitations and recommended applications and uses expressly identified as such shall be binding upon ultrafilter. Samples, descriptions, representations, and other information concerning products in ultrafilter catalogues, advertisements or other promotional materials

or statements of representations made by ultrafilter employees or sales representatives are for general information purpose only and are not binding upon ultrafilter with respect to such warranty provisions.

No employee or sales representative of ultrafilter shall have any authority whatsoever to establish, expand or otherwise modify ultrafilter product specifications, limitations, or recommended applications without ultrafilter giving it's prior express written consent to the customer.

To claim under warranty, the goods must have been installed and continuous maintenance must have taken place in the manner specified in the Operators Handbook. Our service engineers are highly qualified and equipped to assist you in this respect. They are also available to do repairs that may become necessary, in which event they will require an official order before carrying out the work. If such work is to be the subject of a warranty claim, the order should be endorsed 'for consideration under warranty'. Warranty claims must be submitted and shall be processed in accordance with ultrafilter's established warranty claim procedure.

This is the only authorized ultrafilter warranty and is in line with all other expressed or implied warranties or representation including any implied warranties or other obligations on the pad of the manufacturer. The forgoing is the exclusive remedy of any buyer of the manufacturer's product. The maximum damages liability of the manufacturer is the Original purchase price of the product or part. ultrafilter will not be held liable for business interruptions, loss of profits, personal injury, costs of delay or any other special indirect, incidental, special or consequential loss, costs, or damages.